

GT za dermatologiju (2 – 21)

Posude dizajnirane za cuvanje tekućeg dusika za dermatoloske aplikacije

GT opremljen s dva tipa kanistera



Visina kanistera
25 i 50 mm

Kuka fiksira kanister



CRYAL JET maxi

- ▶ U skladu s direktivom European Medical Directive 93/42 EEC
- ▶ Dostupno u 4 velicine: 2, 9, 11 i 21 litar
- ▶ Lagani i jednostavni za upotrebu
- ▶ Kanisteri su zasticeni sigurnosnom kukom
- ▶ Omoguceno rukovanje s obje ruke uz minimalni rizik
- ▶ 6 godina garancija na vakum



CRYAL JET

Cryogenic posude s raspršivacima u spreju za tretiranje tekućim dušikom različiti sprej raspršivaci namijenjeni za tretiranje različitih tipova lezija na koži

Jednostavno punjenje Cryal Jetu potrebom TR pojednostavljenog dekanting sistema **A** TR linija,

Dodaci:



Pojednostavljeni dekanting sistem
Ref.: ACC-DERMATO-1



- ▶ U skladu s direktivom: **European Medical Directive 93/42 EEC**
- ▶ Dostupni u 2 velicine: **0.3 and 0.5 litara**
- ▶ Jednostavni za upotrebu
- ▶ Metalni
- ▶ Odmah raspršuju tekući dušik (sprej)
- ▶ Opremljeni su setom od 6 raspršivaca, dizajniranih za tretiranje točno određenih lezija



Contact probe

Ref.: ACC-CRYALJET-204



Raspršivaci

Ref.: ACC-CRYALJET-105




Cryo rukavice

Ref.: ACC-SECU-15 / 18



Zastitne naočale

Ref.: ACC-SECU-1



Polypro disques

Ref.: ACC-CRYALJET-3

DERMATO

		CRYAL JET				
		GT 2	GT 9	GT 11	GT 21	MAXI
GLAVNE KARAKTERISTIKE:						
Kapacitet	l	2	9.3	12.2	21.5	0.5
Promjer vrata	mm	30	50	50	50	68.5
Staticko vrijeme cuvanja N2	d	25	84	130	225	24
Tezina praznog	kg	1.9	8.2	9.2	13	0.62
Tezina punog	kg	3.5	15.7	19	30.4	0.94
Ukupna tezina	mm	392	450	630	660	285
Promjer kanistera	mm	26	31	31	31	
PRODUCT REFERENCES						
GT DERMATO sa 2 metalna kanistera ¹⁾²⁾		GT2-DERMATO-1	GT9-DERMATO-1	GT11-DERMATO-1	GT21-DERMATO-1	
CRYAL JET MAXI 0.5 l (sa 6raspršivaca)						CRYALJET50-1

1) GT 2 DERMATO:s jednim metalnim kanisterom

2) GT 9 / 11 / 21 DERMATO:sa1, 25 mm visokim kanisterom i 1 50 mm visokim metalnim kanisteromh